## INSTRUCTION FOR USE

Sterile single use administration set/disposable for CO<sub>2</sub> vascular injection.

GB

### INDICATIONS

This product (IS-IN-LL-A) is for use exclusively with CO<sub>2</sub> injector INX-1000A of VIMAP TECHNOLOGIES.

### CONTRAINDICATIONS

- The use of this device and of CO<sub>2</sub> as a contrast medium is contraindicated in the following cases:
  - Supradiaphragmatic arterial procedures (except the case of shunts for haemodialysis and the study of forearm and hands arteries).
  - In case of interatrial or interventricular communication.
  - Pulmonary arteriovenous malformation.
  - Pulmonary insufficiency.
  - Injection in thoracic aorta

### WARNINGS

- Use only this product if the package is unopened, if the aspect of the disposable is clean, not damaged and packaging is sealed and intact.
- Do not reuse or resterilize because you would be exposed to a severe risk of cross contamination or infection because of the ability of some bacteria to form a biofilm on the inner side of the device that could contribute to the failure of decontamination or sterilization process, and a high risk of spore forming by bacteria on a reused product.
- This sterile disposable is for a single use only, and can't be resterilised: It must be replaced every time you perform a procedure on a new patient: risk of contamination of the patient, risk of contamination of the device, risk of failure of one of the components.
- Check to make sure no patient fluids touch the filter. If patient fluids reach the filter, discard the administration set and use a new one.
- Use very cautiously in patients with ischemic bowel or in situations where the gas may cause a "vapor lock"
- Don't use for Supra-diaphragmatic arterial injections, and cerebral administration.
- Relative contraindications to the intravenous use of CO<sub>2</sub> include pulmonary hypertension and chronic obstructive pulmonary disease. As the intravenous injection of CO<sub>2</sub> in diagnostic quantities can increase pulmonary arterial pressure, CO<sub>2</sub> injections must be separated by 3 to 5 minutes of pause to prevent accumulation of CO<sub>2</sub> and subsequent pulmonary artery vapor lock. CO<sub>2</sub> should be used cautiously in patients with patent foramen oval (PFO) or atrial septal defect.
- CO<sub>2</sub> should not be injected into the abdominal aorta in the prone position since the buoyant gas may fill the spinal and lumbar arteries, and cause spinal cord ischemia. Similarly CO<sub>2</sub> should not be injected into the abdominal aorta with the patient's head in an elevated position since the gas can flow in the opposite direction of blood flow, especially in a hypoplastic aorta in children.
- CO<sub>2</sub> should not be injected into the patient undergoing nitrous oxide anesthesia, concurrent use of CO<sub>2</sub> should be avoided since the nitrous oxide can diffuse into the CO<sub>2</sub> bubble, increasing the CO<sub>2</sub> volume significantly. In the venous system, this rapid expansion of CO<sub>2</sub> bubble may result in pulmonary artery vapor lock.
- Each exam has always to be performed under strict monitoring of vital parameters of the patient and more particularly CO<sub>2</sub> concentration in patient's blood and heart rate.

### PRECAUTIONS

- Only connect the IS-IN-LL-A to a needle disposable that is integrated a minimum of one check valve to prevent any accidental contact between the disposable IS-IN-LL-A and patient blood or liquids.
- Pressure, volume and flow must be adequate to the procedure. Please refer to the non-exhaustive list in the INX-1000A CO<sub>2</sub> injector user manual for existing literature, clinical studies or trials.
- An excessive volume of gas in the heart can obstruct the pulmonary outflow tract with severe cardiovascular consequences.


### ALLERGIC REACTIONS

- The VIMAP Administration set is Latex free.
- All plastic parts are disposable, single use parts must not be reused.
- NOTE: Due to anaphylactoid reactions to latex, we strongly recommend that vinyl or non-latex gloves be used during the procedure.

### INSTRUCTIONS FOR USE

- The IS-IN-LL-A disposable is to act as a physical barrier to protect from Virus and Bacteria contamination the device (automatic CO<sub>2</sub> injector) and to protect from this contamination the patient from any virus or bacteria that could be accidentally present in the circuit of the device (between CO<sub>2</sub> supply source and the needle disposable).
- Note 1: The CO<sub>2</sub> supply source (cylinder or wall installation) are not provided by VIMAP TECHNOLOGIES.
- Note 2: The needle disposable connecting the IS-IN-LL-A to the patient is not provided by VIMAP TECHNOLOGIES.

- Please refer to the INX-1000A CO<sub>2</sub> Injector Operator's Manual for additional instructions for use.
- Use pressure adjustment, flow adjustment or injection time to set desired.
- Follow on the injector screen, step by step instructions for connecting and using the consumable.
- For the procedure it is necessary to have this sterile administration set IS-IN-LL-A (provided by VIMAP TECHNOLOGIES), and a Needle disposable with a Check valve/unidirectional output valve (Opening pressure: <=1.5 PSI Back pressure: >=30 PSI), NOT provided by VIMAP TECHNOLOGIES, selected and used by the operator for the injection into the patient artery or vessels.
- As for the Needle disposable to make the injection, the operator must choose the most suitable product available on the market, based on his/her experience and the procedural clinical context.

INX-1000A CO<sub>2</sub> injector  
VIMAP TECHNOLOGIES

IS-IN-LL-A,  
sterile disposable,  
accessory of the INX-1000A.

Provided by  
VIMAP TECHNOLOGIES

Discard the entire administration set after completion of the procedure, in accordance with all Federal, State and Local regulations (These are intended for single use only and sterile).

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