INSTRUCTION FOR USE

Sterile single use administration set/disposable for CO₂ vascular injection.

GB

INDICATIONS

This product (IS-IN-LL-A) is for use exclusively with CO₂ injector INX-1000A of VIMAP TECHNOLOGIES.

CONTRAINDICATIONS

- The use of this device and of CO₂ as a contrast medium is contraindicated in the following cases:
 - Supradiaphragmatic arterial procedures (except the case of shunts for haemodialysis and the study of forearm and hands arteries).
 - In case of interatrial or interventricular communication.
 - Pulmonary arteriovenous malformation.
 - Pulmonary insufficiency.
 - Injection in thoracic aorta

WARNINGS

- Use only this product if the package is unopened, if the aspect of the disposable is clean, not damaged and packaging is sealed and intact.
- Do not reuse or resterilize because you would be exposed to a severe risk of cross contamination or infection because of the ability of some bacteria to form a biofilm on the inner side of the device that could contribute to the failure of decontamination or sterilization process, and a high risk of spore forming by bacteria on a reused product.
- This sterile disposable is for a single use only, and can't be resterilised: It must be replaced every time you perform a procedure on a new patient: risk of contamination of the patient, risk of contamination of the device, risk of failure of one of the components.
- Check to make sure no patient fluids touch the filter. If patient fluids reach the filter, discard the administration set and use a new one.
- Use very cautiously in patients with ischemic bowel or in situations where the gas may cause a "vapor lock"
- Don't use for Supra-diaphragmatic arterial injections, and cerebral administration.
- Relative contraindications to the intravenous use of CO₂ include pulmonary hypertension and chronic obstructive pulmonary disease. As the intravenous injection of CO₂ in diagnostic quantities can increase pulmonary arterial pressure, CO₂ injections must be separated by 3 to 5 minutes of pause to prevent accumulation of CO₂ and subsequent pulmonary artery vapor lock. CO₂ should be used cautiously in patients with patent foramen oval (PFO) or atrial septal defect.
- CO₂ should not be injected into the abdominal aorta in the prone position since the buoyant gas may fill the spinal and lumbar arteries, and cause spinal cord ischemia. Similarly CO₂ should not be injected into the abdominal aorta with the patient's head in an elevated position since the gas can flow in the opposite direction of blood flow, especially in a hypoplastic aorta in children.
- CO₂ should not be injected in the patient undergoing nitrous oxide anesthesia, concurrent use of CO₂ should be avoided since the nitrous oxide can diffuse into the CO₂ bubble, increasing the CO₂ volume significantly. In the venous system, this rapid expansion of CO₂ bubble may result in pulmonary artery vapor lock.
- Each exam has always to be performed under strict monitoring of vital parameters of the patient and more particularly CO₂ concentration in patient's blood and heart rate.

PRECAUTIONS

- Only connect the IS-IN-LL-A to a disposable that is integrated a minimum of one check valve to prevent any accidental contact between the disposable IS-IN-LL-A and patient blood or liquids.
- Pressure, volume and flow must be adequate to the procedure. Please refer to the non-exhaustive list in the INX-1000A CO₂ injector user manual for existing literature, clinical studies or trials.
- An excessive volume of gas in the heart can obstruct the pulmonary outflow tract with severe cardiovascular consequences.

ALLERGIC REACTIONS


- The VIMAP Administration set is Latex free.
- All plastic parts are disposable, single use parts must not be reused.
NOTE: Due to anaphylactoid reactions to latex, we strongly recommend that vinyl or non-latex gloves be used during the procedure.

INSTRUCTIONS FOR USE

- The IS-IN-LL-A disposable is to act as a physical barrier to protect from Virus and Bacteria contamination the device (automatic CO₂ injector) and to protect from this contamination the patient from any virus or bacteria that could be accidentally present in the circuit of the device (between CO₂ supply source and the disposable).
Note 1: The CO₂ supply source (cylinder or wall installation) are not provided by VIMAP TECHNOLOGIES.
Note 2: The disposable connecting the IS-IN-LL-A to the patient is not provided by VIMAP TECHNOLOGIES.

- Follow on the injector screen, step by step instructions for connecting and using the consumable.
- Please refer to the INX-1000A CO₂ Injector Operator's Manual for additional instructions for use.
- Use pressure adjustment, flow adjustment or injection time to set desired.

- For the procedure it is necessary to have this sterile administration set IS-IN-LL-A (provided by VIMAP TECHNOLOGIES), and a disposable with a Check valve/unidirectional output valve (Opening pressure: <=1.5 PSI Back pressure: >=30 PSI), NOT provided by VIMAP TECHNOLOGIES, selected and used by the operator for the injection into the patient artery or vessels.
- As for the disposable to make the injection, the operator must choose the most suitable product available on the market, based on his/her experience and the procedural clinical context.

Discard the entire administration set after completion of the procedure, in accordance with all Federal, State and Local regulations (These are intended for single use only and sterile).

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
IS-IN-LL-A

[GB] **A** Gas input/Connection to the CO₂ injector. **B** Hydrophobic and bacterial filters, Blocks 99.9% of bacteria and viruses, prevents any contamination of the patient or the device. **C** Syringe, to fill with physiological liquid (Physiological liquid not provided by Vimap Technologies). **D** CO₂ sensors / exhaust gas connector, to manage the pressure and to check the concentration of CO₂, inside the disposable. **E** Luer lock compatible Output, to connect the sterile consumable (not provided by VIMAP technologies) that will be in contact with the patient. **F** Clamps, to open or close the tube. **G** One way valve, to prevent any return of patient liquid, and accidental pressure drop of the disposable on the patient side. **H** Stopcock, for directing the fluids according to the steps of the examination. **I** Removable double female luer lock adapter. **J** Plug to remove, before use.

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