## INSTRUCTION FOR USE

Single use non-sterile administration set/disposable for insufflation of the colon during Intussusception Reduction procedure.

<b>INDICATIONS</b>
This product ( <b>AS-3W-H-R35B</b> ), is for use exclusively with CO2 insufflator VMX-1020A of VIMAP TECHNOLOGIES.

- CONTRAINDICATIONS**
- This product should not be used on patients with known or suspected colonic perforation or toxic megacolon. It should not be used within 6 days of large forceps having been used or a “HOT” Biopsy, or snare polypectomy.
  - Do not use this product in a colostomy stoma.
  - Do not use this product following any recent rectal surgery or low rectal anastomosis or when proctitis or other rectal conditions such as inflammatory or neoplastic diseases are suspected.
  - Intussusception Reduction administration sets must not be used for intrauterine insufflation / Distension.

- WARNINGS**
- Use only this product if the package is unopened, if the aspect is clean and packaging is sealed and intact. Do not reuse or sterilize because you would be exposed to a severe risk of cross contamination or infection because of the ability of some bacteria to form a biofilm on the inner side of the device that could contribute to the failure of decontamination or sterilization process, and a high risk of spore forming by bacteria on a reused product.
  - This disposable is for a single use only: It must be replaced every time you perform a procedure on a new patient: risk of contamination of the patient, risk of contamination of the device, risk of failure of one of the components.
  - Care must be taken to avoid over insufflation of the balloon, since with overfilling, or asymmetric filling, displacement of the tip may occur. Such displacement can lead to rectal perforation, or vasovagal reactions, or may cause the balloon to deflate. Over-inflation may cause the inflatable balloon to rupture with possible injury to the patient.
  - The catheter should not be moved unnecessarily once inserted.
  - Forceful or deep insertion may cause tearing or perforation of the rectum.
  - This product is non-sterile, so to avoid any risk of infection, do not use in case of presence of wounds or injuries of colon or patient.

- PRECAUTIONS**
- Inflation of the retention balloon should only be performed when it is deemed necessary and should be performed by a physician, or qualified personnel under the supervision of a physician. Do not inflate the balloon with more than 25 cc water.
  - Do not unnecessarily move the catheter and retention balloon once inserted.
  - The effluent trap should be kept at or below the level of the patient during the procedure. In the unlikely event that it is necessary to elevate the fluid trap above the level of the patient while still connected, the slide clamp must be placed in the completely closed position.
  - Do not use the provided syringe for any other application other than inflating the retention balloon.

- ALLERGIC REACTIONS**
- The VIMAP Administration set is Latex free.
  - All plastic accessories are disposable, single use devices must not be reused or left in the body cavity for an extended period of time. NOTE: Due to anaphylactoid reactions to latex, we strongly recommend that vinyl or non-latex gloves be used during the procedure.

- INSTRUCTIONS FOR USE**
- Please refer to the Colon CO2 Insufflator Operator’s Manual for additional instructions for use.
  - With power turned on, make sure that the CO2 Insufflator is in standby mode and no administration set is connected.
  - Use pressure adjustment to set desired colon pressure.
  - Lubricate the catheter tip and retention Balloon area with water gel lubricant. Do not use Vaseline.
  - Carefully insert the catheter until the proximal end (Base) of the retention balloon is just beyond the anorectal junction. Care must be taken during insertion of the catheter in to the patient not to exert undue pressure to the neuromuscular plexus which can lead to vasovagal reactions and syncopal episodes. Forceful or deep insertion may cause tearing or perforation of the rectum (See WARNINGS).
  - To insufflate the balloon (See PRECAUTIONS), fill the syringe with maximum 25 CC of water, attach the syringe tip to the side port on the catheter and inflate the balloon by discharging the water out of the syringe. Take care to avoid over-inflation of the retention balloon (See WARNINGS). Balloon inflation should only be performed by a physician or qualified personnel under a physician’s supervision. NOTE: Use supplied syringe to inflate the retention balloon only.
  - After inflation, the base of the balloon should be at the level of the anal sphincter for the best retention and maximum safety. If necessary the balloon should be gently pulled back until its proximal end (base) rests on the anal sphincter. NOTE: Position the Administration set so that the tubing is unobstructed.
  - After insertion of the catheter, a discharge of fluid from the colon may be noticed. Assume that the effluent is completely drained, and captured in the effluent trap, before proceeding.
  - Be sure that the fluid trap is not elevated above the level of the patient while the clamp is open (see Precautions)
  - Check to make sure no patient fluids breech the filter.
  - Connect device end of administration set to CO2 Insufflator.
  - Start the flow of CO2, as instructed in the CO2 Insufflator Operator’s Manual.
  - Monitor pressure and volume indicators during the scans.
  - After Completion of the procedure, stop the flow of CO2, as instructed in the CO2 Insufflator Operator’s Manual.
  - Engage Tubing pinch clamp at distal end of administration set.
  - Completely deflate the retention balloon (if inflated) by reattaching the syringe to the side port and pull back the plunger of the syringe. Remove the catheter from the patient.

- Discard the entire administration set after completion of the procedure, in accordance with all Federal, State and Local regulations (These are intended for single use only).

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

# AS-3W-H-R35B

**[GB]** **(A)**: Gas input/connection to the CO<sub>2</sub> insufflator. **(B)**:Hydrophobic and bacterial filters, blocks 99.9% of bacteria and viruses, prevents any contamination of the patient, avoids any contamination of the device in case of reflux. **(C)**: Syringe, use with ambient air, permit to expand the balloon of the catheter. The catheter balloon ensures the tightness of the patient's rectum, and avoids that the probe can come out during the examination. **(D)**: Effluent bag, to store any fluids from the colon. **(E)**: Silicone Foley catheter, in direct contact with the patient, allows to send the gas in the colon, and to collect the possible fluids coming from the colon. **(F)**: Clamps, to open or close the tube.

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