INSTRUCTION FOR USE

Single use non-sterile administration set/disposable for insufflation of the colon during Virtual Colonoscopy (CTC) procedure.

INDICATIONS

This product (**AS-3W-Y-R35A**), is for use exclusively with CO₂ insufflator VMX-1010A de VIMAP TECHNOLOGIES.

CONTRAINDICATIONS

• Do not use in case of : active colonic inflammation, Acute abdominal symptoms, Acute diarrhea or symptoms of acute colitis, Recent acute diverticulitis, Recent colorectal surgery, Symptomatic colon-containing abdominal wall hernia, Recent submucosal endoscopic biopsy or complicated polypectomy/mucosectomy, Known or suspected colonic perforation, Symptomatic or high grade small bowel obstruction, for intrauterine insufflation / Distension.

WARNINGS

- Use only this product if the package is unopened, if the aspect is clean and packaging is sealed and intact. Do not reuse or sterilize because you would be exposed to a severe risk of cross contamination or infection because of the ability of some bacteria to form a biofilm on the inner side of the device that could contribute to the failure of decontamination or sterilization process, and a high risk of spore forming by bacteria on a reused product.
- This disposable is for a single use only: It must be replaced every time you perform a procedure on a new patient: risk of contamination of the patient, risk of contamination of the device, risk of failure of one of the components.
- Care must be taken to avoid over insufflation of the balloon, since with overfilling, or asymmetric filling, displacement of the tip may occur. Such displacement can lead to rectal perforation, or vasovagal reactions, or may cause the balloon to deflate. Over-inflation may cause the inflatable balloon to rupture with possible injury to the patient.
- The catheter should not be moved unnecessarily once inserted.
- Forceful or deep insertion may cause tearing or perforation of the rectum.
- Check to make sure no patient fluids touch the filter. If patient fluids reach the filter, discard the administration set and use a new one.
- This product is non-sterile, so to avoid any risk of infection, do not use in case of presence of wounds or injuries of colon or patient.

PRECAUTIONS

- Inflation of the retention balloon should be performed by a physician, or qualified personnel under the supervision of a physician. Do not inflate the balloon with more than 50 cc water.
- Do not unnecessarily move the catheter and retention balloon once inserted.
- The effluent trap should be kept at or below the level of the patient during the procedure. In the unlikely event that it is necessary to elevate the fluid trap above the level of the patient while still connected, the slide clamp must be placed in the completely closed position.
- Do not use the provided syringe for any other application other than inflating the retention balloon.

ALLERGIC REACTIONS

- The VIMAP Administration set is Latex free.
- All plastic accessoires are disposable, single use devices must not be reused or left in the body cavity for an extended period of time. NOTE: Due to anaphylactoid reactions to latex, we strongly recommend that vinyl or non-latex gloves be used during the procedure.

INSTRUCTIONS FOR USE

- Please refer to the Colon CO₂ Insufflator Operator’s Manual for additional instructions for use.
- Use pressure adjustment to set desired colon pressure.
- Lubricate the catheter tip and retention Balloon area with water gel lubricant. Do not use Vaseline.
- Carefully insert the catheter until the proximal end (Base) of the retention balloon is just beyond the anorectal junction. A «black line» printed on the catheter will indicate this point to the operator: Insert until the «black line» reach patient entrance of the anal canal then you can inflate the balloon. Balloon must always be inflated in the rectal Ampulla area, never in the anal sphincter area.
- To insufflate the balloon (See PRECAUTIONS), fill the syringe with 25-50 CC of water, attach the syringe tip to the side port on the catheter and inflate the balloon by discharging the water out of the syringe. Take care to avoid over-inflation of the retention balloon (See WARNINGS). Balloon inflation should only be performed by a physician or qualified personnel under a physician’s supervision. NOTE: Use supplied syringe to inflate the retention balloon only.
- After inflation, the base of the balloon should be at the level of the anal sphincter for the best retention and maximum safety. If necessary the balloon should be gently pulled back until its proximal end (base) rests on the anal sphincter. NOTE: Position the Administration set so that the tubing is unobstructed, and the effluent trap is below the tubing.
- After insertion of the catheter, a discharge of fluid from the colon may be noticed. Assume that the effluent is completely drained, and captured in the effluent trap, before proceeding.
- Be sure that the fluid trap is not elevated above the level of the patient while the clamp is open (see Precautions)
- Check to make sure no patient fluids breech the filter.
- Connect device end of administration set to CO₂ Insufflator.
- Start the flow of CO₂, as instructed in the CO₂ Insufflator Operator’s Manual.
- Monitor pressure and volume indicators during the scans.
- After Completion of the procedure, stop the flow of CO₂, as instructed in the CO₂ Insufflator Operator’s Manual.
- Completely deflate the retention balloon (if inflated) by reattaching the syringe to the side port and pull back the plunger of the syringe. Remove the catheter from the patient. NOTE: The administration set must be removed from the patient in a timely manner, to allow relief of colonic pressure.

- Discard the entire administration set after completion of the procedure, in accordance with all Federal, State and Local regulations (These are intended for single use only).

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

[GB] **A**: Gas input/connection to the CO₂ insufflator. **B**: Hydrophobic and bacterial filter, blocks 99.9% of bacteria and viruses, prevents any contamination of the patient, avoids any contamination of the device in case of reflux. **C**: Syringe, use with ambient air, permit to expand the balloon of the catheter. The catheter balloon ensures the tightness of the patient’s rectum, and avoids that the probe can come out during the examination. **D**: Effluent bag, to store any fluids from the colon. **E**: Silicone Foley catheter, in direct contact with the patient, allows to send the gas in the colon, and to collect the possible fluids coming from the colon. **F**: Clamps, to open or close the tube. **G**: Blue pant short for the intimacy of the patient. **H**: Release gas tube, to manage the gas over-pressure.

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